



Wearable Breast Pump



Rechargeable & Portable



Low Noise Operation



Anti Backflow and Blockage



Safe For
Baby & Mom



Five Breast Pumping Modes



Massage Mode



Bionic Mode



Breast pumping Mode



Massage+Bionic



Massage+Breast pumping

15 - Gear suction selection in each modes



How to assemble



1. Attach Connector



2. Install Valve



3. Add Diaphragm



3. Add Silicone Shield



4. Attach Milk Collector



5. Snap Into Place



5. Connect Cop

6. Connect Motor Unit

Product Description

1 x Silicone breast shield



1 x **Valve** barbase

1 x Valve base



1 x Silicone Gas-liquid separation



1 x **Main machine**



1 x Duckbill valve



1 x **Pumping bowl**

